




Biocair Packaging Range

Kryotrans Systems available at Biocair (Phase Change Material-based)

System	Temperature range	Specification name	Payload size (cm)	Payload volume	Tare weight	Validation period
KT26 	-15°C to -25°C	KT26-F	25 x 23 x 17.5	10 ℓ	30 kg	approx 48 hrs
	+2°C to +8°C	KT26-RS	25 x 23 x 22	12.6 ℓ	26 kg	
KT70 	-15°C to -25°C	KT70-F	34 x 28.5 x 46	45 ℓ	73 kg	approx 96 hrs
	+2°C to +8°C	KT70-RH	34 x 28.5 x 46	45 ℓ	67 kg	
	+15°C to +25°C	KT70-T	34 x 28.5 x 46	45 ℓ	87 kg	
KT400 	-15°C to -25°C	KT400-F	90 x 50 x 58.5	263 ℓ	205 kg	120 hrs +
	+2°C to +8°C	KT400-RH	90 x 50 x 48.5	218 ℓ	216 kg	
	+15°C to +25°C	KT400-T	90 x 50 x 58.5	263 ℓ	240 kg	



If you need any more information on our shipping systems, please contact your local office.

You can find their numbers on the 'Contact Us' section of our website www.biocair.com.

Or you can email us at enquiries@biocair.com

GLOBAL SERVICE. LOCAL EXPERTISE. PERSONAL ATTENTION.

